Tegaderm™ CHG

Ⓞ	Chlorhexidine Gluconate gel pad device	1
Ⓞ	Tampon de gel au gluconate de chlorhexidine	4
Ⓞ	Chlorhexidingluconat-Gelkissen	7
Ⓞ	Dispositivo con tampone in gel di clorexidina gluconata	10
Ⓞ	Gel con gluconato de clorhexidina	13
Ⓞ	Klorhexidingluconat-gelpude	16
Ⓞ	Gelèpute med klorheksidinglukonat	19
Ⓞ	Klooriheksidiiniglukonaattigeelitynyllinen sidos.....	22
Ⓞ	Penso com dispositivo de camada de gel de gluconato de clorohexidina	25
Ⓞ	Pomôcka s gélovým vankúšikom s obsahom chlórhexidín glukonátu.....	28
Ⓞ	Hlorheksidīna glikonāta gela starplikas izstrādājums.....	31
Ⓞ	Chlorheksidino gliukonato gelio pagalvélės įtaisas.....	34
Ⓞ	Dispozitiv cu pad gel cu gluconat de clorhexidină.....	37
Ⓞ	Засіб із гелевою прокладкою із хлоргексидину глюконату	40
Ⓞ	Jastučić s klorheksidin-gluconat gelom	43
Ⓞ	Подложка с хлоргексидин глюконат гел	46
Ⓞ	Sredstvo sa gel-jastučićem sa hlorheksidin-gluconatom	49
Ⓞ	클로로헥시딘 글루코네이트 젤패드 장치	52

Gluconate gel pad device**Product Description**

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device is used to protect vascular access sites. The Chlorhexidine Gluconate Gel Pad Device is breathable and transparent, allowing continuous site observation.

The transparent film provides an effective barrier against external contamination including fluids (waterproof), bacteria, viruses* and yeast.

The antimicrobial CHG Gel Device consists of an integrated gel pad containing 2% w/w Chlorhexidine Gluconate (CHG), a well-known antiseptic agent with broad spectrum antimicrobial and antifungal activity. The gel absorbs fluid.

**In vitro* testing shows that the transparent film of the gel pad device provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage. The barrier to viruses is due to the physical properties of the dressing, rather than the ancillary properties of CHG.

In vitro testing (time kill and zone of inhibition) demonstrates that the CHG gel pad of the device has an antimicrobial effect against a variety of gram-positive and gram-negative bacteria and yeast.

Intended Use / Intended Purpose:

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device can be used to protect catheter sites. Common applications include protecting intravascular catheters and percutaneous devices. Tegaderm™ CHG is intended to reduce skin colonization and catheter colonization and to suppress regrowth of microorganisms commonly related to bloodstream infections. Tegaderm™ CHG is intended to reduce catheter-related bloodstream infections (CRBSI) in patients with central venous or arterial catheters.

Warnings

DO NOT PUNCTURE CHG GEL WITH NONCORING NEEDLE.

DO NOT USE TEGADERM™ CHG DRESSINGS ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.

THE SAFETY AND EFFECTIVENESS OF TEGADERM™ CHG DRESSINGS HAS NOT BEEN EVALUATED IN CHILDREN UNDER 18 YEARS OF AGE. FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT EARS, EYES, MOUTH OR MUCOUS MEMBRANES.

DO NOT USE THIS PRODUCT ON PATIENTS WITH KNOWN HYPERSENSITIVITY TO CHLORHEXIDINE GLUCONATE. THE USE OF CHLORHEXIDINE GLUCONATE CONTAINING PRODUCTS HAS BEEN REPORTED TO CAUSE IRRITATIONS, SENSITIZATION, AND GENERALIZED ALLERGIC REACTIONS. IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY, AND IF SEVERE, CONTACT A PHYSICIAN.

Hypersensitivity reactions associated with topical use of chlorhexidine gluconate have been reported in several countries. The most serious reactions (including anaphylaxis) have occurred in patients treated with lubricants containing chlorhexidine gluconate, which were used during urinary tract procedures. Caution should be used when using chlorhexidine gluconate containing preparations, and the patient should be observed for the possibility of hypersensitivity reactions.

Precautions

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device should not be placed over infected wounds. It is not intended to be used as a treatment of percutaneous device-related infections. In the case of clinical wound infection, systemic antibacterials should be used if indicated. Any active bleeding at the insertion site should be stabilized before applying the device. Do not stretch the device during application. Mechanical skin trauma may result if the device is applied with tension. The skin should be clean, dry, and free of detergent residue to prevent skin irritation and to ensure good adhesion.

Do not reuse. Reuse may result in compromising product integrity and lead to device failure.

Use care and cover the CHG gel pad device to protect from water when showering.

The CHG gel pad device has not been tested for use during radiation therapy.

Clinical Trial Results: A randomized, controlled clinical trial consisting of 1879 subjects with 4163 central venous and arterial catheter insertion sites was conducted at 11 hospitals (1). Results showed that the use of Tegaderm™ CHG resulted in a statistically significant 60% reduction in the incidence of catheter-related bloodstream infections ($P=0.02$). Study results also demonstrate a statistically significant reduction in skin colonization ($P<0.001$) and catheter colonization ($P<0.0001$) in the chlorhexidine vs. non-chlorhexidine group.

Variable	Non-chlorhexidine vs. chlorhexidine dressings (941 patients/2055 catheters vs. 938 patients/2108 catheters)
Catheter-related bloodstream infection	
Incidence densities (n per 1000 catheter-days)	1.3 vs. 0.5
Hazard ratio	0.402 [0.186 to 0.868], $P=0.02$
Catheter colonization	
Incidence densities (n per 1000 catheter-days)	10.9 vs. 4.3
Hazard ratio	0.412 [0.306 to 0.556], $P<0.0001$


(1) Timsit JF et al (Randomized Controlled Trial of Chlorhexidine Dressing and Highly Adhesive Dressing for Preventing Catheter-Related Infections in Critically Ill Adults *Am. J. Respir. Crit. Care Med.* 2012; 186:1272-1278

Instructions for Use: Refer to packaging figures and those included in package insert. Failure to follow the manufacturer's instructions for use may result in complications including skin irritation and/or maceration.

Site Preparation: Prepare port site according to institution protocol. Clipping of hair at the site may improve CHG gel pad device adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue. Allow all preps and protectants to dry completely before applying CHG gel pad device to prevent skin irritation and to ensure good adhesion.


Any active bleeding at the insertion site should be stabilized before application of CHG gel pad device.

Post Needle Insertion Application of CHG Gel Pad Device (Figures A1 to A5)

1. Open package and remove the sterile CHG gel pad device using aseptic technique.
2. Prepare port site and insert noncoring needle per facility protocol, allowing for a space under the noncoring needle for the CHG gel pad device of about 1/2 cm.
3. Under green tab labeled 1, grasp the soft cloth tabs. Peel and remove top liner and discard (Figure A1).
4. Do not unfold the liner. Pull the soft cloth tabs apart to widen the slit of the gel (Figure A2).
5. Using the soft cloth tabs center the device under the noncoring needle (Figure A3).
6. Lift and grasp the soft cloth tabs and remove liner labeled 2 in direction of arrow (Figure A4).
7. If needed, using the soft cloth tabs, bring slit edges closer together. Do not stretch device during application. Mechanical skin trauma may result if the device is applied with tension (Figure A5).
8. Use firm pressure to smooth down device border and enhance adhesion.
9. See application instructions for 3M™ I.V. Port Dressing.

Simultaneous Application of the CHG Gel Pad Device and Noncoring Needle (Figure B1 to B5)

1. Open package and remove the sterile CHG gel pad device using aseptic technique.
2. Prepare port site and pre-fill noncoring needle system according to facility protocol. Allow all preps and protectants to dry completely before application to prevent skin irritation and to ensure good adhesion.
3. Under green tab labeled 1, grasp the soft cloth tabs. Peel and remove top liner and discard (Figure B1).
4. Do not unfold device liner. Leave protective sheath in place over needle. Orient the CHG gel pad device facing upward and slide needle all the way into base of slit (Figure B2).
5. Wrap the CHG gel pad device around noncoring needle. Remove protective sheath and insert needle into port site per facility protocol (Figure B3).
6. Before removing liner, ensure needle is positioned all the way into base of slit. The CHG gel pad device may be adjusted under noncoring needle by grasping soft cloth tabs. Press gently on soft cloth tabs and remove liner labeled 2 in direction of arrow (Figure B4).
7. If needed, using soft cloth tabs, bring slit edges closer together. Do not stretch device during application. Mechanical skin trauma may result if device is applied with tension (Figure B5).
8. Use firm pressure to smooth down device border and enhance adhesion.
9. See application instructions for 3M™ I.V. Port Dressing.

Site Care:

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the CHG gel pad device, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
2. Inspect the CHG gel pad device daily and change as necessary, in accordance with facility protocol; at a minimum of every 7 days, per Centers for Disease Control and Prevention (CDC) recommendations. CHG gel pad device changes may be needed more frequently with highly exudative sites.
3. Use care and cover the CHG gel pad device to protect from water when showering.

The 3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad Device should be changed as necessary:

- If the device becomes loose, soiled or compromised in any way
- If the site is obscured or no longer visible
- If there is visible drainage outside the CHG gel pad
- If the CHG gel pad device appears to be saturated or overly swollen*

*Note: To test if the CHG gel pad is fully saturated, lightly press down on a corner of the gel pad with your finger. If the CHG gel pad remains displaced once your finger is removed, the CHG device should be changed.

Note: Tegaderm™ CHG gel pad device is not designed to absorb large quantities of blood or fluid.

Removal

CHG Gel Pad Device: See removal instructions for 3M™ I.V. Port Dressing. Remove non-coring needle per facility protocol. Remove CHG gel pad device using soft cloth tabs. A medical adhesive solvent can be used to help remove the device border. If needed, use sterile alcohol swabs or wipes, or sterile solutions (i.e., sterile water or normal saline) to facilitate removal of the gel pad. Continue the low and slow removal method until the CHG gel pad device is completely removed.

Care should be taken not to dislodge catheters or other medical devices when the CHG gel pad device is removed. Support the skin and catheter while removing the CHG gel pad device.

Storage/Shelf Life/Disposal

For best results, store in a cool, dry place. For shelf life, refer to the expiration date on the package.

Sterility of the product is guaranteed unless individual package is damaged or open.

For further information outside of the United States, contact your local 3M representative or contact us at 3M.com and select your country.

How supplied/Ordering Information

Catalog #	CHG Gel Pad Device Size	Average Amount of CHG per CHG Gel Pad Device (mg based on gel pad size)
1665R	CHG Gel Pad Device Dimensions 6.2 cm x 4.9 cm (2 7/16 in x 1 15/16 in)	30

Explanation of Symbols

 Not Made With Natural Rubber Latex

 Caution, see instructions for use

 Do not use if package is damaged or open

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3M and Tegaderm are trademarks of 3M.
© 2020, 3M. All rights reserved.
Made in U.S.A. with
EC REP **3M Deutschland GmbH**