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CE **REF**
Adriatic Health Factory d.o.o.
Nikole Tesle 17
44000 Sisak
CROATIA

Stand der Information: 05/2021 – V7.0

Medizinprodukt: Klasse II b, CE 1250

Instructions for use – ActiMaris® Forte wound irrigation solution

1. Product profile

ActiMaris® Forte wound irrigation solution is made of water with sea salt (3.0 %) and sodium hypochlorite NaOCl (0.2 %). Sodium hypochlorite and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA/VRE tested positive. The oxidative effect of sodium hypochlorite and hypertonic saline are also responsible for the neutralisation of wound odours. The alkaline pH-value of ActiMaris® Forte wound irrigation solution supports due to the provision of OH-ions (as part of a naturally balanced redox system) additionally the release of oxygen. Sea salt contains sodium chloride as its main ingredient (NaCl = 98 %) as well as 2 % other sea salts. The increased salinity of ActiMaris® Forte wound irrigation solution in turn ensures an optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® Forte wound irrigation solution is suitable: for cleaning, humidification, decontamination, biofilm resolution, reduction of swelling and physiological debridement, to create an optimal "micro environment" for:

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, fissures and contusions) and postoperative wounds
- chronic wounds (pressure ulcers, arteriovenous ulcers, diabetic ulcers)
- nerotic, malodorous wounds and ulcerating tumours, even with cavities (hollow spaces)
- thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- venous ulcers

for intra-operative cleaning and irrigation of superficial wounds.

for moistening of bandages and wound dressings.

during dressing changes to loosen encrusted bandages or other encrusted wound dressings.

for the local treatment of the skin and mucous membrane in inflammatory and infectious processes.

for the prevention of infection of the mucous membrane and skin, burns and other wound types.

3. Application instructions

ActiMaris® Forte wound irrigation solution shall be applied at room temperature (not above 30° C). In the event of strong pain, it is suggested that the necessary amount is heated to body temperature in a separate container before use (irrigation solution shall not be heated above 30° C).

Before using ActiMaris® Forte wound irrigation solution, the wounds should first be rinsed and cleaned with ActiMaris® Sensitive wound irrigation solution, in order to remove dirt as well as easily removable coatings. The subsequent deep cleansing and odour reduction will be done with ActiMaris® Forte wound irrigation solution; inert bandages or wound dressings soaked with the wound irrigation solution are used for this purpose (see also Section 9: Interactions).

The application should preferably be in the form of the classic wet/dry phase ("moist balance").

The recommended application time of the wet phase in case of severe fibrinous or purulent, infected or strong malodorous wounds → 5-10 minutes.

The duration and frequency of application should be adapted to the individual wound conditions. The use of the wound irrigation solution is not limited to time and frequency. Necroses and coatings can be removed mechanically, depending on the wound condition. Primary surgical sanitation is advised with thick, adherent, fibrous or necrotic wound coatings.

After the improvement of the wound condition, it is advised to continue further wound cleansing with ActiMaris® Sensitive wound irrigation solution until the wound is visually clean.

ActiMaris® Sensitive wound irrigation solution is recommended for patients with sensitive wounds or skin.

Be careful when applying the wound irrigation solution so that the bottle does not touch the skin or wound of the patient in order to prevent contamination or the spread of micro-organisms.

4. Tissue tolerance and biocompatibility

ActiMaris® Forte wound irrigation solution is suitable for long-term (> 24h) as well as repeated use during bandage changes even with irritated skin and mucous membrane. The application is almost painless and is generally well tolerated by allergy sufferers. Inhibition of granulation and epithelialisation is not to be feared. ActiMaris® Forte wound irrigation solution is dermatologically assessed as safe.

5. Risks following improper use

- Improper use may cause softened wound edges
- Cross-contamination between wounds or patients with incorrect handling of the tube/bottle

6. Side effects

So far, in several years of clinical experience, the following symptoms occasionally occurred:

- Slight temporary burning sensation in case of sensitive wounds

7. Contraindications

Since allergies can never be excluded, ActiMaris® wound irrigation solutions or ActiMaris® wound gel should not be used with known or suspected allergies to any of its ingredients.

8. Restrictions on use

Pregnancy and lactation:

There is no found evidence of genotoxicity, carcinogenicity or embryotoxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast-feeding women, ActiMaris® Forte wound irrigation solution should be used in these cases only after careful medical benefit-risk assessment.

Infants and toddlers:

Due to insufficient experience of long-term exposure, ActiMaris® Forte wound irrigation solution should only be applied after careful medical benefit-risk assessment.

9. Interactions

ActiMaris® Forte wound irrigation solution acts energising and oxidising as a bio-physical compound due to the content of ROS substances (O₂ + OCl⁻), which is why a simultaneous combination with other local, active substances and/or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® Forte wound irrigation solution.

On the other hand, a good compatibility with modern wound dressings such as silver dressings has been shown in clinical trials. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibres, alginates, hydrocolloids etc.

The concomitant use of ActiMaris® Forte wound irrigation solution with ActiMaris® wound gel is possible and useful.

10. General safety instructions

ActiMaris® Forte wound irrigation solution is suitable for external wounds on skin or mucous membrane. Do not use for infusion or injection! Do not swallow!

Use only closed and undamaged containers. Keep out of the reach of children!

It is possible for coloured fabrics to fade after contact with ActiMaris® products. Rinse immediately with water after contact with coloured fabrics.

11. Composition and technical information

Water, sea salt, sodium hypochlorite.

Clear, colourless, alkaline and hypertonic solution with characteristic smell; contains 3.0 % sea salt (Sal Maris) and 0.2 % NaOCl.

12. Storage specification, duration, expiry date

ActiMaris® Forte wound irrigation solution has to be stored at room temperature between +5° and +25° C. Protect from direct sunlight. The shelf life is 30 months. The current expiry date is printed on the packaging. After opening the bottle, the period of use is 3 months; always close the bottle tightly after use.

ActiMaris AG

Manufacturer
ActiMaris AG
Sandgrube 29
9050 Appenzell (SWITZERLAND)
info@actimaris.com
Tel.: +41 71 505 75 25

LOT Batch designation / batch number
REF Order number
Attention, observe instructions for use
Can be used until
Storage instructions (temperature)
Protect from direct sunlight
Manufacturer

CE **REF**
Adriatic Health Factory d.o.o.
Nikole Tesle 17
44000 Sisak
CROATIA

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ActiMaris AG **Fabricant :**