Shiley™

Tracheostomy Tube

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
it	Istruzioni per l'uso
es	Instrucciones de uso
pt	Instruções de uso
nl	Gebruiksaanwijzing
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
el	Οδηγίες χρήσης

Table 1. DCT Dimensions

					
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)	Leak Test Inflation Volume (cc)	Cuff Resting Diameter (mm)
4	5.0	9.4	62	9	20
6	6.4	10.8	74	12	24
8	7.6	12.2	79	15	27
10	8.9	13.8	79	18	29

Table 2. DFEN Dimensions

					
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)	Leak Test Inflation Volume (cc)	Cuff Resting Diameter (mm)
4	5.0	9.4	82	9	20
6	6.4	10.8	74	12	24
8	7.6	12.2	79	15	27
10	8.9	13.8	79	18	29

Table 3. DCFS Dimensions

			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
4	5.0	9.4	62
6	6.4	10.8	74
8	7.6	12.2	79
10	8.9	13.8	79

Table 4. DCFN Dimensions

			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
4	5.0	9.4	62
6	6.4	10.8	74
8	7.6	12.2	79
10	8.9	13.8	79

Table 5.

	DCT	DFEN	DCFS	DCFN
	X	X	X	X
	X	X	X	X
		X	X	X
	X	X	X	X

This instruction insert applies to the following Shiley Tracheostomy products:

DCT	Tracheostomy Tube Cuffed with Disposable Inner Cannula
DFEN	Tracheostomy Tube Cuffed with Disposable Inner Cannula Fenestrated
DCFS	Tracheostomy Tube Cuffless with Disposable Inner Cannula
DCFN	Tracheostomy Tube Cuffless with Disposable Inner Cannula Fenestrated

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

Cuff resting diameter

Tube length

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

Instructions for Use

POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

WARNINGS:

CAREFULLY READ THESE INSTRUCTIONS PRIOR TO USE. THESE INSTRUCTIONS ARE APPLICABLE ONLY TO THE SHILEY DISPOSABLE CANNULA TRACHEOSTOMY TUBES (DCT, DFEN, DCFS, DCFN) AND DISPOSABLE INNER CANNULAE (DIC).

STERILE ONLY IF THE PROTECTIVE BAG IS NOT OPEN, DAMAGED OR BROKEN. DO NOT RESTERILIZE. DO NOT EXPOSE TO TEMPERATURES ABOVE 120°F (49°C).

WHEN USING A LASER IN CONJUNCTION WITH THESE DEVICES, CARE MUST BE TAKEN TO AVOID CONTACT OF THE LASER BEAM WITH THE TUBE. SUCH CONTACT, ESPECIALLY IN THE PRESENCE OF OXYGEN-ENRICHED MIXTURES, COULD RESULT IN RAPID COMBUSTION OF THE TUBE WITH HARMFUL THERMAL EFFECTS AND WITH EMISSION OF CORROSIVE AND TOXIC COMBUSTION PRODUCTS, INCLUDING HYDROCHLORIC ACID (HCl).

DURING AN MRI SCAN THE PILOT BALLOON SHOULD BE SECURED NEAR THE Y CONNECTOR OF THE VENTILATOR CIRCUIT AT LEAST 3 CM FROM THE AREA OF INTEREST TO PREVENT MOVEMENT AND IMAGE DISTORTION.

STERILE. Sterilized with ethylene oxide. For single patient use only.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

CAUTIONS:

The Shiley Tracheostomy Tube and Accessory Products are composed of medical grade plastics which do not contain latex material.

Precautions should be taken when discarding any of the tracheostomy tube devices and accessories. Disposal of medical devices should be made in accordance with applicable national regulations for biologically hazardous waste.

Contents of Packages - Refer to Table 5.

Description - Refer to Tables 1 through 4 for Sizing Dimensions

NOTE: The sizing dimensions as noted on Panel 1 refer to the I.D. as the inside diameter of the inner cannula at the smallest location. O.D. is the outside diameter of the outer cannula (not including cuff features). Length is the distance from the neck plate to the distal tip on the tube's centerline.

The Shiley Disposable Tracheostomy Tubes (DCT, DFEN, DCFS and DCFN) are double cannula tracheostomy tubes with disposable inner cannula and snap-lock connector. The Shiley Tracheostomy Tubes have a radiopaque, biocompatible outer cannula constructed of polyvinyl chloride. A swivel neck plate allows conformity to individual neck anatomies. The Disposable Inner Cannula (DIC) with 15 mm snaplock connector, is translucent for easy inspection and can be used with standard respiratory equipment. The smooth, rounded tip obturator facilitates insertion. The Shiley Disposable Cannula Tracheostomy Tubes (DCT, DFEN, DCFS and DCFN) are available in

four sizes: 4, 6, 8, and 10. Additional Disposable Inner Cannula (DIC) in the same sizes are available separately.

Cuffed models (DCT, DFEN) have a thin wall, high volume, low pressure cuff to minimize pressure in the trachea. When inflated, the cuff conforms to the natural shape of the trachea providing a seal at low intracuff pressure. The cuff inflation line has a luer valve with an integral pilot balloon to indicate cuff inflation.

The red Disposable Decannulation Plug (DDCP) can be used to occlude the proximal end of the outer cannula, forcing the patient to breathe through the fenestrations and upper airway tract during the weaning process. The Disposable Decannulation Plug is available in a universal size that fits all four sizes of the Shiley Disposable Cannula Tracheostomy Tubes (DFEN, DCFS, DCFN) and can be purchased separately.

Indications for Use

The devices are intended for use in providing tracheal access for airway management. The fenestrated devices (DFEN, DCFN) are also indicated when the use of fenestrations is desirable in order to safely and effectively wean a tracheostomy patient. When used in conjunction with the Disposable Decannulation Plug (DDCP), the fenestrated tracheostomy tube can provide a means of weaning and/or phonation for the patient. Use of the Disposable Decannulation Plug with the fenestrated tracheostomy tubes forces air through the fenestrations and around the tube, and into the upper airway and vocal cords.

CAUTION:

The Shiley Tracheostomy Tube is classified as a disposable medical device. Frequent and routine changes of the Tracheostomy Tube are recommended. The manufacturer recommends that the Tracheostomy Tube usage not exceed twenty-nine (29) days. Frequent and routine changes of the Tracheostomy Tube and accessories are recommended and should be evaluated by the attending physician.

The Disposable Inner Cannula (DIC) and the Disposable Decannulation Plug (DDCP) are designed for single use and should not be cleaned or reused.

Instructions

Tube Preparation

1. Selection of trach tube size is left to the discretion of the physician. Patients in the home care environment should be carefully instructed by a home care provider in the proper use and handling of this device.

Pre-Insertion Cuff and Inflation Test

NOTE: Refer to Tables 1 and 2 for Leak Test Inflation Volumes. Inflation volumes are for test purposes only. Consult the physician or Home Health Care Provider for the appropriate inflation volume/pressure when the tube is positioned within the trachea.

2. With the Shiley Disposable Cannula Cuffed models (DCT, DFEN) the cuff and inflation system should be tested for leakage before inserting the tube. This test can be performed as follows: Inflate the cuff with the volume of air indicated in Tables 1 and 2. Then either observe for deflation over several minutes or immerse the tube in sterile saline and observe for air leakage. Deflate the cuff prior to insertion.

Insertion

CAUTION:

To ease insertion and to guard against cuff perforation from sharp edges of cartilage, the cuff should be tapered back. This can be accomplished by first inflating the cuff. Then gently move the cuff away from the distal tip of the outer cannula towards the swivel neck plate as the residual air is removed by deflation. Do not use any sharp instruments such as forceps or hemostats that would damage the cuff when tapering it.

3. Remove the disposable inner cannula and insert the obturator into the outer cannula. The obturator should be fully seated before the tube is inserted into the patient. A thin film of water soluble lubricant can be applied to the outer cannula, cuff, and protruding portion of the obturator to facilitate insertion.
4. After the tracheostomy procedure is performed, insert the tube into the patient's trachea. After the tube is properly in place, remove the obturator immediately.

CAUTIONS:

When using a water soluble lubricant verify that the lubricant does not enter and occlude the tube lumen which may restrict ventilation.

The obturator should be cleaned and then stored in an accessible location near the patient should it be needed for an unscheduled re-intubation.

WARNINGS:

THE SHILEY DISPOSABLE INNER CANNULA PROVIDED IN THIS PACKAGE CAN ONLY BE REPLACED BY A SHILEY DISPOSABLE INNER CANNULA (DIC) OF THE SAME SIZE NUMBER.

A DISPOSABLE INNER CANNULA (DIC) SHOULD BE WITH THE PATIENT AT ALL TIMES. IN THE EVENT OF RESPIRATORY DISTRESS, THE DISPOSABLE INNER CANNULA (DIC) MUST BE INSERTED TO MECHANICALLY VENTILATE THE PATIENT.

DURING AND AFTER ATTACHMENT OF RESPIRATORY OR ANESTHESIA TUBING AND/OR CONNECTORS TO THE DISPOSABLE INNER CANNULA, AVOID APPLICATION OF EXCESSIVE ROTATIONAL, LINEAR, OR ROCKING FORCES ON THE TUBING AND/OR CONNECTORS TO PREVENT ACCIDENTAL DISCONNECTION OF THE DISPOSABLE INNER CANNULA, OR DAMAGE TO THE TRACHEOSTOMY TUBE.

5. Insert the disposable inner cannula with 15 mm connector and lock into position. To lock the disposable inner cannula in place, push connector until both snap-locks

clear the ridged lock on the outer cannula. The outer cannula should always be used with the disposable inner cannula in place unless the Disposable Decannulation Plug (DDCP) is being used.

CAUTION:

To avoid applying pressure against the patient, the tabs on the connector may be held slightly open, and the neck plate may be stabilized with your free hand during the locking procedure.

Verify that the snap-lock connector engages securely. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.

Cuff Inflation

WARNINGS:

UNDER NO CIRCUMSTANCES SHOULD MORE THAN 25 mm OF MERCURY AIR PRESSURE BE USED TO INFLATE THE CUFF. OVER-INFLATION OF THE CUFF MAY CAUSE TRACHEAL DAMAGE AND MAY INHIBIT VENTILATION.

AS A FURTHER PRECAUTION FOR VENTILATOR-DEPENDENT PATIENTS, CUFF INFLATION SHOULD BE CHECKED ON A REGULAR BASIS AND REPLACEMENT TRACHEOSTOMY TUBES SHOULD BE KEPT AT BEDSIDE.

6. Inflate the low pressure cuff by injecting air into the luer valve of the inflation line using a syringe. Selection of a cuff inflation and deflation procedure is left to the discretion of the physician.

CAUTION:

This product is composed of soft materials to conform to tracheal tissue for performance and patient comfort. Simple precautions in handling of the Shiley Cuffed Tracheostomy Tube during insertion and while in place will facilitate proper function and minimize tears and breaks in the inflation system. Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system. It is recommended that the inflation line be maintained in a position allowing for patient mobility without placing tension on the line-to-cannula junction. Prevent lint or other particulates from entering the luer valve of the pilot balloon.

7. Secure the tracheostomy tube to the patient using the neck strap provided.

Cuff Deflation

8. Accumulated secretions above the cuff (DCT, DFEN) may need to be suctioned before deflating the cuff with a syringe, unless suctioning is contraindicated.

9. To deflate the low pressure cuff, withdraw air slowly from the luer valve of the inflation line using a syringe.

Cleaning

WARNINGS:

DO NOT USE SOLUTIONS OR CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED IN THE TABLE BELOW TO CLEAN ANY PART OF THE TRACHEOSTOMY TUBE, AS THIS MAY RESULT IN TUBE DAMAGE.

DO NOT SOAK ANY PART OF THE TUBE IN HYDROGEN PEROXIDE OR ANY OTHER SOLUTION.

<ul style="list-style-type: none">DISPOSABLE INNER CANNULA (DIC)DISPOSABLE DECANNULATION PLUG	DO NOT CLEAN OR REUSE.
<ul style="list-style-type: none">CUFFED OUTER CANNULA	DO NOT CLEAN. RINSE GENTLY WITH STERILE, NORMAL SALINE. AIR DRY.
<ul style="list-style-type: none">CUFFLESS OUTER CANNULANECK PLATEOBTURATOR	<ol style="list-style-type: none">CLEAN WITH HYDROGEN PEROXIDE (HALF-STRENGTH), STERILE NORMAL SALINE, OR WATER AND MILD DETERGENT.AFTER CLEANING THOROUGHLY RINSE WITH STERILE SALINE TO REMOVE ANY CLEANING SOLUTION RESIDUES.

Disposable Decannulation Plug (DDCP)

The red Disposable Decannulation Plug (DDCP) is a universal size that fits all Shiley Disposable models. The disposable decannulation plug occludes the proximal end of the outer cannula, forcing the patient to breathe through the fenestrations and upper airway tract during the weaning process.

WARNINGS:

USE THE DISPOSABLE DECANNULATION PLUG (DDCP) ONLY IN SHILEY DISPOSABLE FENESTRATED LOW PRESSURE CUFFED (DFEN), CUFFLESS (DCFS), AND CUFFLESS FENESTRATED (DCFN) TRACHEOSTOMY TUBES.

PRIOR TO INSERTING THE DISPOSABLE DECANNULATION PLUG (DDCP) MAKE SURE THAT THE FENESTRATIONS ARE NOT OCCLUDED, THAT THE CUFF (DFEN) IS COMPLETELY DEFLATED AND THAT THERE IS SUFFICIENT AIRWAY FOR THE PATIENT.

10. Establish the patency of the patient's upper airway tract. The patient's airway should be cleared by coughing and/or suctioning before plugging the tracheostomy tube with the decannulation plug.
11. Completely deflate the tracheostomy tube cuff (DFEN only).
12. Remove the Disposable Inner Cannula (DIC) with 15mm snap-lock connector.
13. To insert the DDCP, grasp the plug by the pinch clamps of the snap-lock connector and pinch several times to loosen.
14. Holding the pinch clamps, secure one side of the DDCP to the proximal end of the outer cannula.

CAUTION:

To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the locking procedure. Do not attempt to secure both clamps simultaneously.

15. Holding the first side securely, roll the decannulation plug across the connector until the other clamp is in position.
16. Gently but firmly, press the second clamp and let it snap-lock into place. Check for a secure connection.
17. If the patient needs to be mechanically ventilated while the decannulation plug is in use, remove the decannulation plug, insert the Disposable Inner Cannula with 15mm snap-lock connector and then connect it to the ventilator.

WARNING:

IF THE PATIENT HAS DIFFICULTY BREATHING, IMMEDIATELY REMOVE THE DISPOSABLE DECANNULATION PLUG, INSERT THE DISPOSABLE INNER CANNULA AND VERIFY AIRWAY PATENCY.

Removing Tube

18. The disposable inner cannula can be replaced to maintain a patent airway without requiring the entire tube to be changed.
19. Before removing the tracheostomy tube, the cuff (if present) should be completely deflated with a syringe. This will ensure that the cuff passes through the stoma with minimum resistance.

Additional Instructions for Use – Fenestrated Tracheostomy Tubes (DFEN, DCFN)

WARNINGS:

SPECIAL CARE MUST BE EXERCISED TO MAINTAIN FENESTRATION PATENCY, PARTICULARLY FOR PATIENTS WITH THICK, TENACIOUS SECRETIONS THAT MAY BLOCK THE FENESTRATIONS. USE APPROPRIATE HUMIDIFICATION WHEN INDICATED TO ASSIST IN MAINTAINING LUMEN PATENCY.

CARE SHOULD BE TAKEN IN THE SELECTION AND PLACEMENT OF THE DFEN OR DCFN TRACHEOSTOMY TUBE TO MINIMIZE THE COMPLICATIONS THAT MAY ARISE FROM GRANULATION TISSUE. OBSTRUCTION OF THE FENESTRATIONS MAY PREVENT INSERTION OF THE DISPOSABLE INNER CANNULA NEEDED FOR MECHANICAL VENTILATION OR COMPROMISE UPPER AIRWAY ACCESS. IF THIS OCCURS, THE DFEN OR DCFN TRACHEOSTOMY TUBE MAY NEED TO BE REPLACED WITH ANOTHER TYPE OF TUBE.

WHEN USING A FENESTRATED TUBE (DFEN or DCFN) THE DISPOSABLE INNER CANNULA (DIC) MUST BE IN PLACE DURING SUCTIONING TO PREVENT THE SUCTION CATHETER FROM PROTRUDING THROUGH THE FENESTRATIONS OF THE TRACHEOSTOMY TUBE. SUCTIONING WITHOUT THE DISPOSABLE INNER CANNULA (DIC) IN PLACE COULD RESULT IN DAMAGE TO THE TRACHEAL WALL AND COULD CAUSE THE SUCTION CATHETER TO BECOME CAUGHT IN THE FENESTRATIONS.

IF YOU FEEL RESISTANCE WHEN INSERTING THE DISPOSABLE INNER CANNULA (DIC) PAST THE FENESTRATION, DO NOT FORCE THE INNER CANNULA THROUGH THE CANNULA. NOTIFY YOUR HEALTH CARE PROVIDER IMMEDIATELY.

NEVER USE THE FENESTRATED TRACHEOSTOMY TUBES FOR MECHANICAL VENTILATION UNLESS THE DISPOSABLE INNER CANNULA WITH 15mm CONNECTOR IS IN PLACE TO PREVENT LEAKAGE THROUGH THE FENESTRATIONS.

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
EC|REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.