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INSTRUCTION FOR USE

NOVO CM 20 Smoothbore S


Catheter Mount

REF N22222020

NOVO

IFU_Rev. 01_202209

Bestandteile Components

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EN | INSTRUCTION FOR USE

Airway connector for conduction of respiratory gases between the breathing system and the patient interface (e.g. endotracheal tube, mask) of a mechanically ventilated patient. For single use only in children, adolescents and adults.

1. NOTICE ON REPORTING REQUIREMENTS

Please inform the manufacturer and the competent authority of the Member State in the event of any serious incident that has occurred in relation to the device.

2. Application by healthcare professionals or users trained by healthcare professionals.

CONTRAINDICATIONS

The device is used only as directed. Improper use of it will harm the patient.

WARNING NOTICES

- The medical device must not be reused or reprocessed.
- Modifying the medical device may cause damage and endanger the patient.
- The gas sampling port must only be used for monitoring breathing gas and anesthesia gases.
- The port for suction or bronchoscopy should only be used for this function. Do not introduce gases or liquids through it.
- Regularly check the inside of the medical device for liquid or secretion accumulation.
- If includes: a Luer Lock connector must only be used for gas monitoring.
- The device is sterilized with Ethylene Oxide. Sterile conditions are guaranteed only if the package is not damaged.

Novo Klinik-Service GmbH
Technologiepark West – Haus 15 |
Zum Frenser Feld 1
50127 Bergheim | Germany

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