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 Hangzhou Shanyou Medical Equipment Co., Ltd.
No. 138, Louta Development Zone, Guancun Village,
Louta Town, Xiaoshan District, Hangzhou, Zhejiang
311266, China

 Importeur
Novo Klinik-Service GmbH
Technologiepark West - Haus 15
Zum Frenser Feld 1
50127 Bergheim | Germany
Tel. +49 2271 99404-0
Fax +49 2271 99404-100
info@its-novo.de | its-novo.de

  Prolinx GmbH
Brehmstr. 56 | 40239 Düsseldorf, Germany
Tel: +49-2131-4051968-0 | E-mail: info@eulinx.de

Endotracheal Tube Standard | with cuff

REF

- IC-30mC IC-35mC IC-40mC IC-45mC IC-50mC IC-55mC
 IC-60mC IC-65mC IC-70mC IC-75mC IC-80mC IC-85mC

Quantity:
1 pc/pack
20 pcs/box

PRODUCT DESCRIPTION

The Endotracheal tube is inserted into the trachea through glottis to establish and maintain an artificial airway during anesthesia and emergency treatment.

INTENDED USERS

Nursing staff, medical professionals

INDICATIONS

The product is intended for use in the emergency department, ICU, or other departments for anesthesia, first aid or resuscitation etc.

CONTRAINDIKATIONS

- The product should be contraindicated in patients allergic to the product material.
- The product should be used cautiously in patients with laryngeal edema, acute airway inflammation and throat abscess.
- The product should be used cautiously in patients with thoracic aortic aneurysm compression of trachea and severe bleeding constitution. The Endotracheal tube is inserted into the trachea through glottis to establish and maintain an artificial airway during anesthesia and emergency treatment.

APPLICATION DESCRIPTION

- Choose the right tube size according to the patient's age, gender, height, weight etc.
- Using aseptic technique, test the cuff by fully inflating it and noting that it does not leak.
- Prior to intubation, deflate the cuff completely.
- After intubation, take out the stylet and inflate the cuff using the minimum volume of air required to provide an effective seal.
- Immediately after cuff inflation, auscultate both lung fields. If breath sounds are diminished over one lung field or absent over one or both fields, adjust the tube as required.
- Endotracheal tube placement should be confirmed by viewing the position of the tube tip with a chest radiograph.

WARNING NOTICES

- Be sure all connections are secure.
- Keep the tracheal tube unobstructed and timely aspirate secretions.
- The patients with tracheal intubation should receive oral cleaning and care regularly.
- Strengthen the management of tube airway's temperature and humidity.
- During the indwelling time, release and re-inflate the cuff regularly to alleviate the pressure of trachea for a long time.
- Endotracheal intubation is intended for use in a short time. If further treatment is needed, other forms of intubation is recommended.

LOT NUMBER

The LOT number is printed on the packaging.

EXPIRY DATE

The expiry date is imprinted on the product and outer packaging.

STORAGE CONDITIONS

Do not store in direct sunlight, at extreme temperatures or humidity. Store in dry conditions, in a non-corrosive atmosphere and in a well-ventilated room.

HANDLING CONDITIONS

In case of damaged packaging or inappropriate storage the product SHOULD NOT be used and has to be destroyed.

STERILE CONDITION

The product is sterilised by ethylenoxide. It is delivered in a single sterile barrier system (○) with a protective outer packaging (⊗).

SINGLE USE NOTICE RESP. NOTICE ON REUSE

The product is a disposable product. Do not reuse.

REPROCESSING NOTICE

The product SHOULD NOT be reprocessed.

DISPOSAL

The product can be disposed with normal clinical waste or according to local regulations. Special guidelines need to be taken into consideration.

NOTICE ON REPORTING REQUIREMENTS

Please inform the manufacturer, the distributor and the importer in all cases of serious incidents in connection with the product.

Hangzhou Shanyou Medical Equipment Co., Ltd.
No. 138, Louta Development Zone, Guancun Village,
Louta Town, Xiaoshan District, Hangzhou, Zhejiang
311266, China

Prolinx GmbH
Brehmstr. 56 | 40239 Düsseldorf, Germany
Tel: +49-2131-4051968-0 | E-mail: info@eulinx.de

Importer
Novo Klinik-Service GmbH
Technologiepark West - Haus 15
Zum Frenser Feld 1
50127 Bergheim | Germany
Tel. +49 2271 99404-0
Fax +49 2271 99404-100
info@its-novo.de | its-novo.de