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Novo SealPad | Adhesive Tracheostomy wound dressing

1. PRODUCT DESCRIPTION

The Novo SealPad Wound Dressing is a conformable, sterile, self adhesive wound dressing composed of (i) a silicone wound contact layer (0.2mm thick) and (ii) a polyurethane foam pad and (iii) a water-proof pu film backing and (iv) a wound location device for placement of a patient's or carer's finger.

Novo SealPad absorbs exudate from the wound and is impermeable to bacteria. The absorption process is characterized by swelling of the polyurethane foam pad in contact with wound exudates through the perforated silicone wound contact layer, providing moist conditions conducive to natural healing. The soft silicone wound contact layer is non-adherent to the wound so that newly formed tissue is not traumatised when Novo SealPad is changed or removed.

2. INDICATIONS

Novo SealPad may be used for the management of the neck wound after removal of a tracheostomy tube, until healing is complete. Post decannulation site of tracheostomy.

3. CONTRAINDICATIONS | SAFETY INFORMATION

Novo SealPad is a sterile product and should not be used if the pouch is damaged or has previously been opened. The device cannot be resterilised.

Novo SealPad should not be used on:

ulcers caused by chronic infections (tuberculosis, deep mycotic infections, syphilis)

arteriopathy stage IV, bites and 3rd degree burns.

Novo SealPad should not be used when there is a known sensitivity to silicone.

Novo SealPad should not be used together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.

Wounds showing clinical signs of infection (temperature, pus, inflammatory signs) should be treated under medical control before use of Novo SealPad can be resumed.

4. INFORMATION FOR USE

Due to the nature of this dressing it is suitable for the closing of a tracheostomy site, as it will withstand the pressure from airflow attempting to pass through the remaining stoma, it can manage the presence of the existing secretions and is also easily removed for necessary changes without damaging surrounding skin or the wound edges.

Novo SealPad is very simple to apply, requiring no special skills or equipment. Dressing changes should normally occur when the polyurethane foam pad has been fully-absorbed; therefore, the interval between dressing changes will depend entirely upon the state of the wound. In the early post decannulation period, and on large or moderately exuding wounds, changes may be required more frequently but may be reduced as healing progresses and exudate levels are reduced (see Section 4.3).

4.1 Skin Preparation

- Cleanse the wound with sterile saline and sterile swabs.
- Dry the skin surrounding the wound.
- Remove any remains of creams or greasy substances.

4.2 Dressing Application

- Select the appropriate Novo SealPad size that will completely cover the wound surface, ensuring a 2 to 3 cm margin beyond the edges of the wound.
- Remove the protective liner from the Novo SealPad and align the central button directly/centrally over the stoma site. Having aligned the Novo SealPad apply to the skin, gently pressing the edges of the dressing. No additional fixation tape is necessary.
- The alignment of the button over the stoma site allows more accurate application of finger pressure by the patient/carer
- Ensure the patient/carer is aware of the clinical requirements to apply pressure to the button as per instructions in section 5.

4.3 Dressing Changes

Novo SealPad is easy to remove due to the soft silicone layer when a whitish blister has formed, and approaches the edges of the dressing. This makes for more comfortable dressing changes with little or no associated trauma, helping to produce improved wound healing.

Novo SealPad should be changed as required, depending on the amount of exudate from the wound. Where loss of adhesion around the wound occurs, the dressing should be changed immediately. Where leakage occurs, the dressings should be changed immediately.

- Gently remove Novo SealPad.
- Remove any remaining debris or gel from the wound, using sterile compresses.
- Follow procedure 4.1 a) to 4.2 c) to apply new dressing.

Novo SealPad is a single use device and once removed **MUST NOT** be reused or reapplied due to the risk of cross infection or re-infection.

Dispose to clinical waste or if at home to domestic waste.

5. SPECIAL NOTES & CAUTIONS

Novo SealPad is a dynamic dressing - in that it requires patient interaction- the patient must be informed of the importance when using Novo SealPad of firm finger pressure to the Central button whenever they speak or cough. Failure to achieve this patient compliance significantly reduces the benefit of this unique dressing.

In the unlikely instance that there is some tissue necrosis around the tracheostomy wound then the wound may initially appear to increase in size in the early stages of Novo SealPad treatment. This is normal and occurs as any wound debris is removed from the edges of the wound.

In the management of moderately to heavily exuding wounds, Novo SealPad can only make the overlying environment more conducive to healing. There are cases where healing is impaired as a result of underlying conditions; in these instances, Novo SealPad alone may make little or no progress, and suitable treatment of the underlying conditions will be necessary as well. Therefore, if after 4-6 weeks of Novo SealPad treatment, there has been no improvement then, in line with accepted wound management practice, the original diagnosis and overall therapy should be reassessed.

Novo SealPad should be left in place as long as possible, but not more than 5 days in order to prevent trauma to the fragile newly formed tissue and to reduce cross contamination through frequent dressing changes.

Thick necroses should be removed before applying Novo SealPad.

Novo SealPad Wound Dressings are intended for single use only.

6. STORAGE | STERILISATION | SYMBOLS ON LABELLING:

Store Novo SealPad dressings away from direct sunlight at ambient temperature (5-40°C) and humidity.

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CE 0123

7. PRESENTATIONS:

Novo SealPad is available in the following sizes:

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Huizhou Foryou Medical Devices Co., Ltd.
CORP. GmbH (Europe)
Novo Klinik-Service GmbH