
## TuBo-Clip / TuBo-Clip XL

REF  
20-01  
20-05  
  
Article  
TuBo-Clip  
TuBo-Clip XL

Patientenzielgruppe: Intubierte Patienten  
Patient target group: Intubated Patients

**Novo Klinik-Service GmbH**  
Technologiepark West - Haus 15  
Zum Freienseer Feld 11 50127 Bergheim | Germany  
Tel. +49 2271 994-04-0 | Fax +49 2271 994-04-100  
info@novo-med.de | novo-med.de

IFU\_Rev. 04\_202211

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### EN | Instructions for use

One-piece, latex-free tube holder for securing endotracheal tubes (oral/nasal), double lumen tubes, laryngeal tubes and laryngeal masks. It is applied in adolescents and adults.

### 1. Useful life

Change in case of contamination / defect / according to hygiene guidelines. Maximum operating life 24 hours.

### 2. Notice on reporting requirements

Please inform the manufacturer and the competent authority of the Member State in the event of any serious incident that has occurred in relation to the device.

### 3. Notice on medical devices

Medical device of Class I. Accessory product for endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks.

### 4. Contraindications

Do not use in case of injuries or accesses in the operating area of the tube holder or allergies to the materials used.

### 5. Warning notices

- The tube holder is applied on the connected tube or on the connected mask.
- Apply the tube holder in two and constantly ensure a secure fit of the tube or mask.
- The safe hold of the fixation has to be checked regularly.
- If the fixation is too loose, extubation can be caused. A slight tension on the tube holder is optimal.
- Make sure that the patient's skin is not exposed to the hook and loop strip of the tube holder, as this may cause irritations and open wounds.
- When the patient is moving, the fixation around the back of the head may slip and become loose thereby. This can result in extubation.
- The product is a disposable product. The reuse of disposables can result in infection and/or effects on the functioning.

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