
REF	Article
20-10	TuBo <sup>®</sup> -Fix
20-10K	TuBo <sup>®</sup> -Fix 1
20-10-02	TuBo <sup>®</sup> -Fix 2
20-10-03	TuBo <sup>®</sup> -Fix 3
20-10-05	TuBo <sup>®</sup> -Fix 4
20-10-TS	TuBo <sup>®</sup> -Fix 5
20-11	TuBo <sup>®</sup> -Fix Mittelteil
20-12	TuBo <sup>®</sup> -Fix Mittelteil XL
20-11-03	TuBo <sup>®</sup> -Fix Mittelteil 3
20-12-03	TuBo <sup>®</sup> -Fix Mittelteil 3 XL
20-11-05	TuBo <sup>®</sup> -Fix Mittelteil 5
20-12-05	TuBo <sup>®</sup> -Fix 5 Mittelteil 5 XL
20-15	TuBo <sup>®</sup> -Fix Nackenpolster
20-16	TuBo <sup>®</sup> -Fix Nackenpolster XL

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## EN | Instructions for use

Latex-free tube holder for securing oral endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks. It is applied in adolescents and adults.

1. Useful life  
Change in case of contamination / defect / according to hygienic guidelines. Maximum operating life 24 hours.
2. Notice on reporting requirements  
Please inform the manufacturer and the competent authority of the Member State in the event of any serious incident that has occurred in relation to the device.
3. Notice on medical devices  
Medical device of Class I. Accessory product for endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks.
4. Contraindications  
Do not use in case of injuries or accesss in the operating area of the tube holder or allergies to the materials used.
5. Warning notices
  - The tube holder is applied on the connected tube or on the connected mask.
  - Apply the tube holder in twos and constantly ensure a secure fit of the tube or mask.
  - The safe hold of the fixation has to be checked regularly.
  - If the fixation is too loose, extubation can be caused. A slight tension on the tube holder is optimal.
  - Make sure that the patient's skin is not exposed to the hook and loop strip of the tube holder, as this may cause irritations and open wounds.
  - Make sure that the adhesive tape does not come into contact with the patient's skin.
  - The product is a disposable product. The reuse of disposables can result in infection and/or effects on the functioning.

Application by healthcare professionals or users trained by healthcare professionals. Decision on use is made by the treating doctor.

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